

RNA extraction from frogs

Preparation the day before the extraction

Materials needed

- Beads
- 2 ml tubes for tissue
- 100% & 70% ethanol
- (2) β -mercaptoethanol⁽¹⁾ and fume hood access
- Qiagen RNA mini kit
- Filter tips (1000 μ l)⁽²⁾
- Glass beaker for β -mercaptoethanol waste (labelled)
- Glass beaker for other waste (labelled)

Preparation

Once per kit: Add 44 μ l ethanol to RPE buffer, and tick the box on the lid.

Collect tissue in RNAlater.

Extraction day

Setup⁽³⁾

1. Clean the bench, trays, pipetes, centrifuge and gloves with 70% ethanol.
2. Prepare 2 ml tubes with sample names.
3. Prepare 1 long tray with pink RNEasy mini spin columns, label with numbers.
4. Prepare tubes with 600 μ l 70% **ethanol**.
5. Prepare RLT buffer in 2 ml tubes: per 1 sample 600 μ l **RLT** + 6 μ l **β -mercaptoethanol**. Work in fume hood.
6. Prepare glass beaker in fume hood (RLT buffer + β -mercaptoethanol waste) and glass beaker for other waste⁽⁴⁾.

⁽¹⁾ In tissue culture room 4°C.

⁽²⁾ In qPCR room, stock in end of same corridor.

⁽³⁾ Since RNA is sensitive, the more time is spent preparing the lab before the RNA is in a sensitive stage, the better the extraction quality.

⁽⁴⁾ Could wash down sink but they react badly with bleach, so better to evaporate in fume hood.

Extraction

Work in fume hood until step 10.

1. Place tissue in labelled 2 ml tubes. 30 mg max, about 3 mm³ cube. Add 600 μ l prepared **RLT** and 1-2 glass **beads**. Use a light-coloured tray.
2. Take photo of ordered tubes in tray. Write the names in order in the lab book.
3. Disrupt tissue in beadbeater, 2 min, 25 Hz.
4. Centrifuge for 3 min at full speed.
5. Use filter tips. Remove supernatant to tubes with 70% **ethanol** and mix immediately by pipetting. Transfer 700 μ l to pink RNEasy mini spin columns.
6. Centrifuge at 8000 RCF for 15 s.
7. Discard flow through in fume hood beaker for β -mercaptoethanol waste⁽⁵⁾. Transfer remaining supernatant to same spin columns.
8. Centrifuge at 8,000 RCF for 15 s.
9. Discard flow through in fume hood beaker for β -mercaptoethanol waste.
10. Add 700 μ l **RW1** solution⁽⁶⁾.
11. Spin for 15 sec at 8,000 RCF.
12. Discard flow-through in other waste glass beaker. Add 500 μ l **RPE**.
13. Spin for 15 sec at 8,000 RCF.
14. Discard flow-through in other waste glass beaker. Add 500 μ l **RPE**.
15. Spin for 2 min at 8,000 RCF.
16. Prepare 1 tray with RNase free kit tubes for RNA elution. Label with names from photo you took.
17. Discard flow-through in other waste glass beaker. Transfer column to RNase-free tubes from kit. Add 50 μ l **RNase-free H₂O** from kit, directly on the membrane.

⁽⁵⁾ Leave beaker to evaporate in fume hood. Wash salt remains.

⁽⁶⁾ For all these washes, close the lid and roll and invert the tubes to achieve a thorough wash. DNase treatment done here if needed.

18. Clean centrifuge lid and tube-holding area with 70% ethanol.
19. Spin for 1 min at 8,000 RCF⁽⁷⁾.
20. (Optional) Pipette 3 μ l to additional tube "Q" to be used for nanodrop, qubit and tape-station. Keep on ice. Order in a box and transfer samples to -80°C until use.
21. If not already done so, update lab book with extraction names.
22. Store samples at -80°C.

⁽⁷⁾ Single elution yields 200 - 800 ng/ μ l in 50 μ l from g28 embryos.